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The active terminal units AGSS type 1L are the connection points between the receiving and disposal systems of the anaesthetic gases where the operator makes the connections and disconnections by using a suitable gas-specific probes. The vacuum is generated by an exhaust ejector Venturi located inside the terminal unit and driven by compressed air.

This feature allows a great versatility of the terminal unit, as it's enough to drive it with compressed air to obtain the requested vacuum, without using any external vacuum sources, such property makes it suitable for small and big health structures.

The terminal unit is composed of two parts: the base block, intended to be connected to the disposal system by means of welding nut; and the socket, connected to the base by a threaded connection and that contains the connection point where the medical equipments are linked to the gas specific probes.

The terminal units respect the requirements of the harmonised standard ISO 9170-2. In order to comply also the requirements of the ISO 7396-2 standard, concerning the anaesthetic gas disposal system, this device includes a correct work visual green indicator of the exhaust ejector Venturi, pneumatically driven.

The device is preset in order to meet the flow rate requirements of the ISO 7396-2 standard, it is also possible by a specific integrated device, to set a different flow-rate.

This terminal unit can be installed on bedhead units or supplied with wall embedded or external kit for wall mounting.

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Outlet and inlet for anaesthetic gas scavenging system – venturi model

Technical Data

INLET NOMINAL PRESSURE: 4,5 bar
MIN INLET NOMINAL PRESSURE: 3 bar
MAX INLET NOMINAL PRESSURE: 5 bar
DEPRESSION GENERATED WITH AIR INLET PRESSURE AT 5 BAR: -0,025 bar
FLOW-RATE WITH AIR INLET PRESSURE AT 5 BAR: 80 l/min
MAX AIR CONSUMPTION (with 50 l/min outlet flow-rate): 30 l/min
OPERATING TEMPERATURE: +10 °C / +40 °C
STORAGE TEMPERATURE: -20 °C / +60 °C
WEIGHT: 0,90 Kg
AVAILABLE BOX:
 - wall-embedded kit
 - external kit
STANDARD: ISO 9170-2 Tipo 1L; ISO 7396-2

Product Codes

AGSS outlet – Venturi model

9320-2028V AGSS outlet – Venturi model

AGSS probe

8634-2029 AGSS outlet probe

Boxes


9322-1600 Wall-embedded outlet kit

9318-1600 External outlet kit

LEGEND

- | | |
|---|------------------------|
| 1 | Socket |
| 2 | Base block |
| 3 | 4 mm Polyurethane hose |
| 4 | Operating lamp |
| 5 | Non return valve |
| 6 | Flow-rate regulation |
| 7 | Plastic hose |
| 8 | Hose connection |


AGSS outlet drawing- Venturi model


Outlet and inlet for anaesthetic gas scavenging system – **venturi model**

Page 2 of 2



AGSS external outlet kit – Venturi model

AGSS wall-embedded outlet kit – Venturi model

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