

## Specification

Model	VS-10H/VS-10M/VS-10S
Dimensions	191(H) x 92 (W) x 112 (D) mm
Weight	About 0.25kg (including the battery)
Power Adapter	Input voltage: AC 100-240V 50/60Hz Input power: 25VA. DC output voltage: 5V 2A
Battery	Built-in battery: 3.6V 3400mAh Continuous operation duration of the lithium battery: 4 hours Lithium battery charge time: not longer than 4 hours
Rotation Angle of Display Screen	Maximum vertical rotation angle: 140°±10° Maximum horizontal rotation angle: 270°±10°
Screen	3.5 inches touchscreen with resolution of 640x960 pixels
Water Proof Level	IP66
Photo Record	Photo can be recorded in JPEG format, up to 10,000 photos can be recorded
Video Record	Video can be recorded in MPEG-4 format, up to 2 hours video can be recorded
Media Storage	Not less than 5GB
Data Port	USB Type C, HDMI*, and Wi-Fi** are available
Anti-Fog	Immediate anti-fog technology is supported with no need for heating of the anti-fog coating on the laryngoscope blade
Blades	Single use blade
Frequency for WiFi**	2.4GHz

Note: \*For VS-10H & VS-10M, \*\*For VS-10H only



## Smaller Size, Bigger Vision

VS-10 Series Video Laryngoscope

# VS-10 Series Video Laryngoscope



## Vision with Safety

- Smart Power Management  
Remaining Time Display
- IP66, Waterproof Design
- Instant anti-fog technology  
No heating time required

## Vision with Ease

- 3.5"HD touchscreen: a step closer to reality
- Horizontal rotation angle: 0°~270°
- Vertical rotation angle: 0°~140°

## Vision For Everyone

All 5 blades share a general handle  
Disposable blades offer a better infection control

