


300 is widely used in chest, Abdominal surgery, ENT, gynaecology and obstetrics, urology and orthopedics etc.

Hydraulic lift by foot pedal, head operated movements.

The base and column cover all made by premium 304 stainless steel.

Tabletop is made of composite laminate for X-ray, makes a high definition image.



Model	Site	Length	Width	Tabletop lowest height	Elevation range
300		≥2020mm	500mm	≤760mm	250mm

**PARAMETER** 300 Manual Hydraulic Operating Table

Movements		Model: 300
Tabletop trendelenburg	Reverse Trendelenburg	≥25°
	Trendelenburg	≥25°
Lateral Tilt (Left And Right)		≥20°
Head plate	Up	≥45°
	Down	≥70°
Back plate	Up	≥75°
	Down	≥15°
Leg plate	Up	≥15°
	Down	≥90°
	Outward	≥90°
Kidney Bridge Elevation		≥110 mm