

COMEN

COMEN Share with the World



C11

Fetal & Maternal Monitor



Shenzhen Comen Medical Instruments Co., Ltd.

P/N: EN-C11-4P-20200302-V1.2

Add: No.2 of FIYTA Timepiece Building, Nanhuan Avenue, Matian Sub-district, Guangming District, Shenzhen, 518106, Guangdong, China.
Tel:+86-755-2640 8879 400-700-9488 Fax:+86-755-2643 1232 Website: en.comen.com E-mail: info@szcomen.com

CE 1639

C11 Fetal & Maternal Monitor



- Folding-up screen, easy for observation.



- 5.6" LED screen.
- Built-in high-capacity Li-ion battery, above 2 hours working time.
- 150mm wide printer. Special extractable paper slot, easy for paper replacement and observation.
- Support 3cm/Min, 2cm/Min, 1cm/Min
- real-time printing.
- Support timing printing.
- 3-level audio/visual alarms.



Support data storage and software upgrading by external USB disk. External SD card optional, for data transferring.



Portable hidden handle.

- Advanced patient files management system.
- Manual/automatic FM counting, automatic fetal status scoring system.
- FHR signal strength indicator, twins cross channel verification.
- 60 hours real-time waveform review
- Timing monitoring function, avoid overtime monitoring.
- Wire & wireless central monitoring station connection.

