



300C

Electric Operating Table

<p>Optional 1006 Orthopedic Traction</p>	<p>Optional 1006 Tibia Bone-Nailing Position</p>	<p>Optional Lateral Arm Support</p>
<p>Optional Waste Basin for Gynecology</p>	<p>Anorectum Operation Position</p>	<p>E.N.T. Operation Position</p>
<p>Built-in Kindey Bridge</p>	<p>Tabletop Sliding</p>	
<p>Optional Head Rest</p>	<p>Elevation Up/Down</p>	
<p>Lateral Tilt Left/Right</p>	<p>Trendelenburg/Reverse Trendelenburg</p>	

Item	Model	ET300C
Tabletop Length		2070mm
Tabletop Width		550mm
Tabletop Minimum Height		700mm
Tabletop Maximum Height		1000mm
Trendelenburg/Reverse Trendelenburg		25°/25°
Lateral Tilt		15°/15°
Head Rest(Up/Down)		45°/90°
Foot Plate(Up/Down)		15°/90°
Foot Plate Horizontal Positioning		90°
Back Plate(up/down)		75°/20°
Kidney Bridge Elevation		110mm
Horizontal Sliding		300mm
Voltage		220V/110V
Frenquency		50Hz/60Hz
Power Capacity		1.0Kw
Battery		Yes
Mattress		Memory Mattress
Accessory		PU By Moulding
Main Material		304 Stainless Steel
Maximum Load Capacity		250KG
Warranty		2 Years

Movements Diagram of ET300C

