



Model	Site	Length	Width	Tabletop lowest height	Elevation range
MT600		≥2000mm	500mm	≤670mm	300mm

**PARAMETER** MT600 Manual Hydraulic Operating Tables

Movements		Model:MT600
Tabletop trendelenburg	Reverse Trendelenburg	≥25°
	Trendelenburg	≥30°
Lateral Tilt (Left And Right)		≥20°
Head plate:	Up	≥60°
	Down	≥75°
Back plate:	Up	≥70°
	Down	≥30°
Leg plate:	Up	≥15°
	Down	≥90°
	Outward	≥90°
Kidney Bridge Elevation		≥ 120 mm